

Starkey EDGE AI

Uncompromising sound. Unrivaled ease.



ITER mRICR RIC RT **RIC 312** (NW) (NW)



Lasting Power, **Lifelong Clarity**

- Longest-lasting rechargeable RIC (up to 51 hours)
- Onboard battery holds three charges
- 10 minutes of charging delivers approximately three hours of use

Tech Tiers Channels Bands	Premium 24	Advanced 20	Select 16
Edge Mode+			
User Control Personalization App Control	•	••	0 0
Speech Optimization			
Consonant Enhancement Voice Al	•	0	0
Sound Manager			
Ear-to-Ear Technology Quiet	•	•	•
Machine Noise Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Speech in Noise Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Transient Noise Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Wind Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Directional Processing*			
Adaptive Dynamic	•	•	•



StarLink Wireless Accessories



Mini Remote Microphone



Remote Microphone +



Remote Control 2.0



Table Microphone



TV Streamer

Compatible Chargers



StarLink Charger 2.0



StarLink Premium Charger 2.0**



StarLink Custom Charger 2.0



StarLink Premium Custom Charger 2.0**

**Premium: Wireless charging case with three additional full charges

Key Features	24 Channels	20 Channels	16 Channels
Voice Al	•		
Translate	•		
Multiflex Tinnitus Pro	•	•	
Transcribe	•	•	
Edge Mode+	•	•	•
Smart Assistant	•	•	•
2-Way Audio & Find My Phone	•	•	•
Hearing Health & Activity Tracking	•	•	•
Fall Alert	•	•	•
Reminders & Self Check	•	•	•
Auto Sleep	•	•	•
Tap Control	•	•	•
SnapFit 2.0 & Custom Cased SnapFit 2.0 Receivers	•	•	•
TeleHear full remote fit & REM Target Match	•	•	•
StarLink wireless accessory compatibility	•	•	•
My Starkey & Hear Share app compatibility	•	•	•





Na CI Salt Mist Test






Extended
Reliability Testing**



Salt Brine Immersion Test**





Pro8 HydraShield™





